


Anti-Carbonic Anhydrase 2 (Ca2) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P00918

Item Number: CM0012142



Introduction

Carbonic Anhydrase 2 (CA2) Rabbit Monoclonal Antibody for WB, IHC-P, ELISA. Reacts with Human, Mouse, Rat. Targets a synthetic peptide within aa 161-260 of CA2 (P00918). Ideal for pH regulation and bone resorption research.

[Learn More](#)

Parameter	Value
Product Name	Carbonic Anhydrase 2 (CA2) Rabbit mAb
Remarks/Alias	CAC; CAII; Car2; CA-II; HEL-76; HEL-S-282; Carbonic Anhydrase 2 (CA2)
Species	Human
GeneID	760
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 161-260 of human Carbonic Anhydrase 2 (CA2) (P00918)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00918
Protein Weight	29kDa
Shipping	Ice bag