

Anti-Cardiac Troponin T Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P45379

Item Number: CM0012844



Introduction

Rabbit polyclonal antibody against human cardiac troponin T (TNNT2), validated for WB, IF-P, IHC-P, ELISA; cross-reacts with mouse and rat; protein weight 31-36 kDa.

[Learn More](#)

Parameter	Value
Product Name	Cardiac troponin T (TNNT2) Rabbit pAb
Remarks/Alias	CMH2; RCM3; TnTC; cTnT; CMD1D; CMPD2; LVNC6; Cardiac troponin T (TNNT2)
Species	Human
Gene ID (Human)	7139
Gene ID	7139
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Cardiac troponin T (Cardiac troponin T (TNNT2)) (NP_001001431.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P45379
Protein Weight	31-36 kDa
Shipping	Ice bag