


Anti-Cas9 Rabbit Monoclonal Antibody For Wb, Elisa - Q99Zw2

Item Number: CM0004241



Introduction

Rabbit monoclonal antibody against Cas9 (SpCas9) validated for WB and ELISA. Immunogen spans aa 450-730 of *S. pyogenes* Cas9. Predicted cross-reactivity with human and other species. Ideal for CRISPR research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Cas9 Rabbit mAb
Remarks/Alias	CAS9; Cas9
Species	Staphylococcus aureus
GeneID (Human)	CAS9
GeneID	901176
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 450-730 of streptococcus pyogenes m1 Cas9. (NP_269215.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Other (Wide Range Predicted)
SWISS	Q99ZW2
Protein Weight	158kDa
Shipping	Ice bag