

Anti-Casp8Ap2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9UKL3

Item Number: CM0017406



Introduction

Rabbit polyclonal antibody against human CASP8AP2 (FLASH), validated for WB and ELISA. Reacts with human. Targets amino acids 1733–1982. Ideal for studying TNF-induced NF- κ B activation, apoptosis, and histone mRNA processing.

[Learn More](#)

Parameter	Value
Product Name	CASP8AP2 Rabbit pAb
Remarks/Alias	CED-4; FLASH; RIP25; CASP8AP2
Species	Human
GeneID (human)	9994
GeneID	9994
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1733-1982 of human CASP8AP2 (NP_036247.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UKL3
Protein Weight	223kDa
Shipping	Ice bag