

Anti-Caspase-6 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P55212

Item Number: CM0002635

Introduction

Anti-Caspase-6 rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Detects human and mouse caspase-6 (33 kDa), a protease mediating apoptosis, PANoptosis, and antiviral defense. Research-use IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	Caspase-6 Rabbit mAb
Remarks/Alias	MCH2; CSP-6; caspase-6; Caspase-6
Species	Human
Gene ID (Human)	839
Gene ID	839
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human Caspase-6 (P55212).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P55212
Protein Weight	33kDa
Shipping	Ice bag