

Anti-Caspase-8 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - O89110

Item Number: CM0009034

Introduction

Rabbit monoclonal anti-Caspase-8 antibody for mouse, validated in WB, IP, IF/ICC, and ELISA. Targets key apoptosis initiator (UniProt O89110) for cell death and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Caspase-8 Rabbit mAb
Remarks/Alias	MACH; Mch5; FLICE; CASP-8; Caspase-8
Species	Mouse
GeneID (Human)	841
GeneID	12370
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of mouse Caspase 8 (NP_033942.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	O89110
Protein Weight	55 kDa
Shipping	Ice bag