

Anti-Caspase-9 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P55211

Item Number: CM0013241



Introduction

Caspase-9 Rabbit Polyclonal Antibody (pAb) is validated for Western blot, immunohistochemistry, immunofluorescence, and ELISA. It recognizes human, mouse, and rat caspase-9, making it a valuable tool for apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	Caspase-9 Rabbit pAb
Remarks/Alias	MCH6, APAF3, APAF-3, PPP1R56, ICE-LAP6, Caspase-9
Species	Human
Gene ID (Human)	842
Gene ID	842
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 139-416 of Caspase-9 (NP_001220.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55211
Protein Weight	46kDa
Shipping	Ice bag