

Anti-Casr Polyclonal Antibody For Wb, Ihc-P, Elisa - P41180

Item Number: CM0018461



Introduction

CASR Rabbit Polyclonal Antibody (PAb) validated for WB, IHC-P, ELISA. Detects human, mouse, rat CASR. Targets human CASR (P41180). Suitable for calcium homeostasis, kidney, and parathyroid research.

[Learn More](#)

Parameter	Value
Product Name	CASR Rabbit pAb
Remarks/Alias	CAR; FHH; FIH; HHC; EIG8; HHC1; NSHPT; PCAR1; hCasR; GPRC2A; HYPOC1; CASR
Species	Human
GeneID (Human)	846
GeneID	846
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-536 of human CASR (NP_000379.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41180
Protein Weight	121kDa
Shipping	Ice bag