


Anti-Cathepsin E (Ctse) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P14091

Item Number: CM0025433



Introduction

Cathepsin E (CTSE) rabbit polyclonal antibody validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying immune function, antigen processing, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin E (CTSE) Rabbit pAb
Remarks/Alias	CATE; Cathepsin E (CTSE)
Species	Human
Gene ID (Human)	1510
Gene ID	1510
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 132-244 of human Cathepsin E (CTSE) (P14091)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14091
Protein Weight	43kDa
Shipping	Ice bag