

Anti-Cathepsin L Monoclonal Antibody For Wb, Ihc-P, Elisa - P07711

Item Number: CM0006462



Introduction

Cathepsin L Rabbit mAb for WB, IHC-P, ELISA. Recognizes human Cathepsin L (38kDa). Important thiol protease in lysosomes, antigen processing, and viral entry. Research use only.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin L Rabbit mAb
Remarks/Alias	MEP; CATL; CTSL1; Cathepsin L
Species	Human
GeneID	1514
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human Cathepsin L (NP_001903.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P07711
Protein Weight	38kDa
Shipping	Ice bag