

Anti-Ccl4 Rabbit Polyclonal Antibody For Wb, Elisa - P14097

Item Number: CM0021492

Introduction

CCL4 (MIP-1-beta) rabbit polyclonal antibody validated for WB and ELISA in mouse. Recognizes a 10 kDa secreted chemokine involved in inflammation and chemotaxis. Ideal for immunodetection studies.

[Learn More](#)

Parameter	Value
Product Name	CCL4 Rabbit pAb
Remarks/Alias	Act-2; Mip1b; Scya4; MIP-1B; AT744.1
Species	Mouse
GeneID (Human)	6351
GeneID	20303
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-92 of mouse CCL4 (NP_038680.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P14097
Protein Weight	10kDa
Shipping	Ice bag