

Anti-Ccr5 Polyclonal Antibody (Rabbit) For If-P And Elisa - P51681

Item Number: CM0022605



Introduction

High-quality rabbit polyclonal antibody against human CCR5, suitable for immunofluorescence and ELISA. Detects mouse CCR5. Supports research on inflammation, immune cell migration, and HIV entry.

[Learn More](#)

Parameter	Specification
Product Name	CCR5 Rabbit pAb
Remarks/Alias	CKR5; CCR-5; CD195; CKR-5; CCCKR5; CMKBR5; IDDM22; CC-CKR-5; CCR5
Species	Human
Gene ID (Human)	1234
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CCR5 (NP_000570.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse
SWISS	P51681
Protein Weight	41kDa
Shipping	Ice bag