

Anti-Ccr6 Polyclonal Antibody For Wb And Elisa - P51684

Item Number: CM0022718



Introduction

High-quality rabbit polyclonal antibody against human CCR6 (C-C chemokine receptor type 6). Validated for Western blot and ELISA. Ideal for investigating immune cell chemotaxis, inflammation, and mucosal immunity.

[Learn More](#)

Parameter	Value
Product Name	CCR6 Rabbit pAb
Remarks/Alias	BN-1; DCR2; DRY6; CCR-6; CD196; CKRL3; GPR29; CKR-L3; CMKBR6; GPRCY4; STRL22; CC-CKR-6; C-C CKR-6; CCR6
Species	Human
GeneID (Human)	1235
GeneID	1235
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human CCR6 (NP_004358.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P51684
Protein Weight	42kDa
Shipping	Ice bag