

Anti-Cct6B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92526

Item Number: CM0014398

Introduction

High-quality rabbit polyclonal antibody targeting human CCT6B, validated for Western Blot, Immunofluorescence, and ELISA. Detects endogenous levels in human, mouse, and rat samples. Ideal for chaperonin complex and testis-specific protein research.

[Learn More](#)

Specification	Value
Product Name	CCT6B Rabbit pAb
Remarks/Alias	Cctz2; CCTZ-2; TSA303; CCT-zeta-2; TCP-1-zeta-2; CCT6B
Species	Human
Gene ID (Human)	10693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-340 of human CCT6B (NP_006575.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92526
Protein Weight	58kDa
Shipping	Ice bag