

Anti-Cd103 Rabbit Polyclonal Antibody For Wb, Elisa - Q60677

Item Number: CM0023887

Introduction

Rabbit polyclonal antibody targeting mouse CD103 (integrin alpha-E, alpha M290). Suitable for Western blot and ELISA applications. Reacts with mouse. Useful for cell adhesion and mucosal immunity research.

[Learn More](#)

Parameter	Value
Product Name	CD103 Rabbit pAb
Remarks/Alias	CD103; aM290; alpha-E1; A530055J10; alpha-M290
Species	Mouse
Gene ID	16407
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 191-385 of mouse CD103 (NP_032425.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q60677
Protein Weight	129kDa
Shipping	Ice bag