

Anti-Cd105/Endoglin Monoclonal Antibody For Elisa, Wb, Fc - P17813

Item Number: CM0011081



Introduction

Rabbit monoclonal antibody against human CD105/Endoglin, validated for ELISA, WB, FC. Targets angiogenesis regulator; reacts specifically with human endoglin, a TGF-beta coreceptor. Ideal for vascular biology and IVD research.

[Learn More](#)

Parameter	Value
Product Name	CD105/Endoglin Rabbit mAb
Remarks/Alias	END; HHT1; ORW1
Species	Human
GeneID	2022
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-337 of human CD105/Endoglin (NP_001108225.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	ELISA, WB, FC
Cross Reactivity	Human
SWISS	P17813
Protein Weight	67kDa/70kDa
Shipping	Ice bag