

Anti-Cd10/Mme Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P08473

Item Number: CM0007559



Introduction

Rabbit monoclonal antibody against human CD10/MME (Neprilysin, CALLA). Validated for WB, IHC-P, IF-P, ELISA. Reacts with human, mouse, rat. Suitable for studying neprilysin function in peptide degradation, cardiovascular regulation, and disease.

[Learn More](#)

Parameter	Value
Product Name	CD10/MME Rabbit mAb
Remarks/Alias	NEP; SFE; CD10; CALLA; CMT2T; SCA43; CD10/MME
Species	Human
Gene ID (Human)	4311
Gene ID	4311
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of human CD10/MME (NP_000893.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08473
Protein Weight	86kDa
Shipping	Ice bag