

Anti-Cd117/C-Kit Rabbit Monoclonal Antibody For Wb, Elisa - P10721

Item Number: CM0006639



Introduction

High-quality anti-CD117/c-Kit rabbit monoclonal antibody validated for WB, ELISA. Targets human KIT receptor (P10721), key in cell survival, hematopoiesis, stem cell maintenance, and melanogenesis. Ideal for cancer and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	CD117/c-Kit Rabbit mAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC; CD117/c-Kit
Species	Human
GeneID (Human)	3815
GeneID	3815
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-520 of human c-Kit (NP_000213.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag