


Anti-Cd123/IL-3R Alpha Polyclonal Antibody (Rabbit) For Wb And Elisa - P26951

Item Number: CM0019944



Introduction

Rabbit polyclonal antibody targeting human CD123/IL-3RA, the ligand-binding alpha subunit of the IL-3 receptor. Validated for WB and ELISA, shows cross-reactivity with mouse. Ideal for hematopoiesis, leukemia stem cell research, and signal transduction studies.

[Learn More](#)

Parameter	Value
Product Name	CD123 Rabbit pAb
Remarks/Alias	IL3R, CD123, IL3RX, IL3RY, IL3RAY, hIL-3Ra, IL-3R-alpha
Species	Human
GeneID	3563
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-305 of human CD123 (NP_002174.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P26951
Protein Weight	43kDa
Shipping	Ice bag