

Anti-Cd138/Syndecan-1 Rabbit Monoclonal Antibody For Ihc-P, Elisa - P18827

Item Number: CM0007267



Introduction

This rabbit monoclonal antibody detects human CD138/Syndecan-1, a cell surface proteoglycan. Suitable for IHC-P and ELISA applications in research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	CD138/Syndecan-1 Rabbit mAb
Remarks/Alias	SDC; CD138; SYND1; syndecan; CD138/Syndecan-1
Species	Human
GeneID (human)	6382
Immunogen	Recombinant protein of human CD138/Syndecan-1
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P18827
Protein Weight	32kDa
Shipping	Ice bag