

Anti-Cd138/Syndecan-1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P18827

Item Number: CM0017893



Introduction

High-quality rabbit polyclonal antibody against human CD138/Syndecan-1 (UniProt P18827). Validated for WB, IHC-P, ELISA; cross-reacts with human and mouse. Ideal for exosome and matrix research.

[Learn More](#)

Parameter	Value
Product Name	CD138/Syndecan-1 Rabbit pAb
Remarks/Alias	SDC; CD138; SYND1; syndecan; CD138/Syndecan-1
Species	Human
Gene ID (Human)	6382
Gene ID	6382
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-240 of human CD138/Syndecan-1 (NP_002988.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P18827
Protein Weight	32kDa
Shipping	Ice bag