

Anti-Cd144 Rabbit Polyclonal Antibody For Wb, Elisa - P55284

Item Number: CM0023665



Introduction

Rabbit polyclonal antibody targeting mouse CD144 (VE-cadherin), raised against recombinant protein. Validated for WB and ELISA applications. Ideal for endothelial cell adhesion and junction studies.

[Learn More](#)

Parameter	Value
Product Name	CD144 Rabbit pAb
Remarks/Alias	7B4; Vec; VECD; Cd144; VEcad; VE-Cad; CD144
Species	Mouse
GeneID (Human)	1003
GeneID	12562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 46-592 of mouse CD144 (NP_033998.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P55284
Protein Weight	88kDa
Shipping	Ice bag