

Anti-Cd152/Ctla-4 Polyclonal Antibody For Wb, Ihc-P, Elisa - P16410

Item Number: CM0022324



Introduction

Detect human CTLA-4 (CD152) with this rabbit polyclonal antibody, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Recombinant immunogen targeting the extracellular domain, ideal for immune checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	CD152/CTLA-4 Rabbit pAb
Remarks/Alias	CD; GSE; GRD4; ALP55; CD152; CTLA-4; IDDM12; CELIAC3; CD152/CTLA-4
Species	Human
GeneID (Human)	1493
GeneID	1493
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16410
Protein Weight	25kDa
Shipping	Ice bag