


Anti-Cd155/Pvr Rabbit Monoclonal Antibody For Wb, Elisa - P15151

Item Number: CM0010490



Introduction

Rabbit monoclonal antibody against human CD155/PVR (Poliovirus receptor), validated for Western blot and ELISA. Suitable for research on NK cell biology, tumor immunology, and viral infection mechanisms.

[Learn More](#)

Parameter	Value
Product Name	CD155/PVR Rabbit mAb
Remarks/Alias	PVS; HVED; CD155; NECL5; TAGE4; Necl-5; CD155/PVR
Species	Human
Gene ID (Human)	5817
Gene ID	5817
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CD155/PVR (P15151).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15151
Protein Weight	45kDa
Shipping	Ice bag