

Anti-Cd16 Rabbit Monoclonal Antibody For Wb, Elisa - P08637 / O75015

Item Number: CM0005111



Introduction

A rabbit monoclonal antibody recognizing human CD16 (FcγRIIIA), with cross-reactivity to mouse, validated for WB and ELISA applications, ideal for ADCC and immune function studies.

[Learn More](#)

Parameter	Value
Product Name	CD16 Rabbit mAb
Remarks/Alias	CD16; FCG3; CD16A; FCGR3; IGFR3; IMD20; FCR-10; FCRIII; CD16-II; FCGRIII; FCRIIIA; FcGRIIIA
Species	Human
Gene ID (Human)	2214/2215
Gene ID	2214/2215
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human CD16. (P08637).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08637/O75015
Protein Weight	29KDa/26KDa
Shipping	Ice bag