

Anti-Cd168/Rhamm Rabbit Polyclonal Antibody (Ko Validated) For Wb, If, Elisa - O75330

Item Number: CM0017041

Introduction

Rabbit polyclonal antibody targeting human CD168/RHAMM, KO validated. Recombinant immunogen (aa1-300). Detects the hyaluronan receptor in WB, IF/ICC, and ELISA. Supports research on cell motility, metastasis, and signaling.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CD168/RHAMM Rabbit pAb
Remarks/Alias	CD168; IHABP; RHAMM; MM
Species	Human
Gene ID (Human)	3161
Gene ID	3161
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CD168/RHAMM (NP_001136028.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O75330
Protein Weight	84kDa
Shipping	Ice bag