

Anti-Cd177 Polyclonal Antibody For Wb, Elisa - Q8N6Q3

Item Number: CM0026565

Introduction

CD177 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting the CD177 antigen (NB1/PRV1) expressed on human neutrophils. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CD177 Rabbit pAb
Remarks/Alias	NB1; PRV1; HNA2A; PRV-1; HNA-2a; NB1 GP; CD177
Species	Human
Gene ID (Human)	57126
Gene ID	57126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-285 of human CD177 (NP_065139.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N6Q3
Protein Weight	51kDa
Shipping	Ice bag