


Anti-Cd2 Monoclonal Antibody For Wb, Fc, Elisa - P06729

Item Number: CM0004026



Introduction

This CD2 rabbit monoclonal antibody recognizes human CD2, a T-cell surface antigen involved in adhesion and signaling. Suitable for WB, FC, and ELISA applications. Protein weight ~39kDa. Ideal for immunological research.

[Learn More](#)

Parameter	Value
Product Name	CD2 Rabbit mAb
Remarks/Alias	T11; SRBC; LFA-2
Species	Human
GeneID (Human)	914
GeneID	914
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-209 of human CD2 (NP_001758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag