

Anti-Cd200 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P41217

Item Number: CM0002815

Introduction

Rabbit monoclonal antibody suitable for WB, IHC-P, and ELISA targeting human CD200 (OX-2 membrane glycoprotein). Supports immunology research on immune regulation and T-cell costimulation. Swiss-Prot: P41217.

[Learn More](#)

Parameter	Value
Product Name	CD200 Rabbit mAb
Remarks/Alias	MRC; MOX1; MOX2; OX-2; CD200
Species	Human
Gene ID (Human)	4345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-232 of human CD200 (NP_001352780.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P41217
Protein Weight	31kDa
Shipping	Ice bag