

Anti-Cd248/Tem1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hcu0

Item Number: CM0015372

Introduction

Rabbit polyclonal antibody targeting human CD248/TEM1, validated for WB and ELISA. Recognizes the tumor endothelial marker endosialin involved in angiogenesis and immune modulation. Suitable for research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	CD248 Rabbit pAb
Remarks/Alias	TEM1; CD164L1; CD248
Species	Human
GeneID (Human)	57124
GeneID	57124
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-150 of human CD248 (NP_065137.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HCU0
Protein Weight	81kDa
Shipping	Ice bag