


## Anti-Cd304/Neuropilin-1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P97333

Item Number: CM0004214



### Introduction

Rabbit monoclonal antibody targeting mouse CD304/Neuropilin-1, validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for cardiovascular development, angiogenesis, and neuroscience research.

[Learn More](#)

| Parameter        | Specification                                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD304/Neuropilin-1 Rabbit mAb                                                                                                   |
| Remarks/Alias    | Nrp; NP-1; Npn1; NPN-1; C530029I03                                                                                              |
| Species          | Mouse                                                                                                                           |
| Gene ID (Human)  | 8829                                                                                                                            |
| Gene ID          | 18186                                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-856 of mouse CD304/Neuropilin-1 (NP_032763.2). |
| Source           | Rabbit                                                                                                                          |
| Category         | Monoclonal Antibodies                                                                                                           |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                               |
| SWISS            | P97333                                                                                                                          |
| Protein Weight   | 103kDa                                                                                                                          |
| Shipping         | Ice bag                                                                                                                         |