


Anti-Cd324/E-Cadherin Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P09803

Item Number: CM0009780



Introduction

High-quality rabbit monoclonal antibody targeting CD324/E-Cadherin (UM, UVO, Ecad, ARC-1, E-cad, L-CAM). Suitable for WB, IF-P, IHC-P, ELISA in mouse and rat samples. Predicted molecular weight 98kDa. UniProt: P09803.

[Learn More](#)

Parameter	Value
Product Name	CD324/E-Cadherin Rabbit PolymAb®
Remarks/Alias	Um; UVO; Ecad; ARC-1; E-cad; L-CAM
Species	Mouse
GeneID (Human)	999
GeneID	12550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 157-369 of mouse CD324/E-Cadherin (NP_033994.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P09803
Protein Weight	98kDa
Shipping	Ice bag