

Anti-Cd326/Epcam Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q99Jw5

Item Number: CM0005819

Introduction

Rabbit monoclonal antibody targeting mouse CD326/EpCAM, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with human and rat. Suitable for studying epithelial adhesion, mucosal immunity, and tumor biology.

[Learn More](#)

Parameter	Value
Product Name	CD326/EpCAM Rabbit mAb
Remarks/Alias	EGP; Ly74; gp40; CD326; EGP-2; TROP1; Egp314; Ep-CAM; EpCAM1; Tacsd1; GA733-2; Tacstd1
Species	Mouse
GeneID (Human)	4072
GeneID	17075
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-266 of mouse CD326/EpCAM (NP_032558.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99JW5
Protein Weight	35 kDa
Shipping	Ice bag