


Anti-Cd328 Polyclonal Antibody For Wb, Elisa - Q9Y286

Item Number: CM0021593



Introduction

High-quality rabbit polyclonal antibody against human CD328 (SIGLEC7), validated for WB and ELISA. Ideal for natural killer cell inhibitory receptor studies. Recombinant immunogen; cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CD328 Rabbit pAb
Remarks/Alias	p75; QA79; AIRM1; CD328; AIRM-1; CDw328; D-siglec; SIGLEC-7; SIGLECP2; SIGLEC19P; p75/AIRM1
Species	Human
Gene ID (Human)	27036
Gene ID	27036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-144 of human CD328. (NP_055200.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y286
Protein Weight	51kDa
Shipping	Ice bag