


Anti-Cd352/Slamf6 Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - Q96Du3

Item Number: CM0006337



Introduction

Rabbit monoclonal antibody targeting human CD352/SLAMF6 (SLAM family member 6), suitable for WB, IHC-P, flow cytometry, and ELISA. Detects 24 and 37 kDa isoforms. Useful for studying NK cell activation and T-cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	CD352/SLAMF6 Rabbit mAb
Remarks/Alias	SLAMF6; CD352; KALI; KALib; Ly108; NTB-A; NTBA; SF2000; SLAM family member 6; CD352/SLAMF6
Species	Human
GeneID (Human)	114836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-215 of human CD352/SLAMF6 (NP_001171643.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	Q96DU3
Protein Weight	24 kDa/37 kDa
Shipping	Ice bag