

Anti-Cd38 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28907

Item Number: CM0013586



Introduction

Rabbit polyclonal antibody against human CD38, validated for WB, IHC-P, IF/ICC, ELISA. Detects endogenous CD38 in human, mouse, and rat samples. Ideal for studying ADP-ribosyl cyclase activity, calcium signaling, and immune cell function.

[Learn More](#)

Parameter	Value
Product Name	CD38 Rabbit pAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1; CD38
Species	Human
GeneID (Human)	952
GeneID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (NP_001766.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28907
Protein Weight	34kDa
Shipping	Ice bag