

Anti-Cd3E Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - P07766

Item Number: CM0002554

Introduction

Rabbit monoclonal antibody against human CD3E, validated for WB, IF, IHC, and ELISA. Part of TCR-CD3 complex essential for T-cell activation. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	CD3E Rabbit mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon; CD3E
Species	Human
Gene ID	916
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, mIHC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag