


Anti-Cd40 Rabbit Monoclonal Antibody For Wb And Elisa - P25942

Item Number: CM0006296



Introduction

Rabbit monoclonal antibody against human CD40 (TNFRSF5). Validated for WB and ELISA. Detects B-cell surface antigen CD40, a key immune receptor. Ideal for cancer and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD40 Rabbit mAb
Remarks/Alias	p50; Bp50; CDW40; TNFRSF5; CD40
Species	Human
GeneID (Human)	958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-193 of human CD40. (NP_001241.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P25942
Protein Weight	31kDa
Shipping	Ice bag