

Anti-Cd45 Rabbit Monoclonal Antibody For Fc, Elisa - P08575

Item Number: CM0012581



Introduction

CD45 rabbit monoclonal antibody (clone mAb) specifically recognizes human CD45 (PTPRC, leukocyte common antigen). Validated for flow cytometry and ELISA. Suitable for immunological studies and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD45 Rabbit mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	P08575
Protein Weight	131kda/136kDa/141kDa/143kDa/147kDa
Shipping	Ice bag