


Anti-Cd45 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P08575

Item Number: CM0005809



Introduction

CD45 (PTPRC) rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Targets human receptor-type tyrosine-protein phosphatase C, essential for T-cell activation and leukocyte function. Suitable for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CD45 Rabbit mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 194-390 of human CD45 (NP_002829.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08575
Protein Weight	147kDa
Shipping	Ice bag