

Anti-Cd62P/P-Selectin Monoclonal Antibody For Wb, Ihc-P, Elisa - P16109

Item Number: CM0005786



Introduction

High-quality anti-CD62P/P-selectin rabbit monoclonal antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Specific to the 91kDa protein, ideal for inflammation and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	CD62P/P-selectin Rabbit mAb
Remarks/Alias	CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM; CD62P/P-selectin
Species	Human
GeneID (Human)	6403
GeneID	6403
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 731-830 of human CD62P/P-selectin (P16109)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16109
Protein Weight	91kDa
Shipping	Ice bag