

Anti-Cd68 Rabbit Monoclonal Antibody For Wb, If-P, Fc (Intra), Elisa - P34810

Item Number: CM0003685

Introduction

Recombinant rabbit monoclonal antibody to human CD68 (macrosialin, GP110). Validated for WB, IF-P, FC (intra), and ELISA. Reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	CD68 Rabbit mAb
Remarks/Alias	GP110; LAMP4; SCARD1; CD68
Species	Human
GeneID (Human)	968
GeneID	968
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD68 (NP_001242.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34810
Protein Weight	37kDa
Shipping	Ice bag