

Anti-Cd68 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P34810

Item Number: CM0015335

Introduction

Rabbit polyclonal antibody against human CD68 (macrosialin, GP110, LAMP4). Validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for macrophage biology and endolysosomal studies.

[Learn More](#)

Parameter	Value
Product Name	CD68 Rabbit pAb
Remarks/Alias	GP110; LAMP4; SCARD1; CD68
Species	Human
Gene ID (Human)	968
Gene ID	968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-200 of human CD68 (NP_001242.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34810
Protein Weight	37kDa
Shipping	Ice bag