

Anti-Cd69 Rabbit Monoclonal Antibody For Wb, Elisa - Q07108

Item Number: CM0006932



Introduction

Recombinant rabbit monoclonal antibody against human CD69 (SWISS: Q07108) for western blot and ELISA. Detects CD69 in human samples and cross-reacts with rat. Immunogen corresponds to amino acids 62-199.

[Learn More](#)

Parameter	Value
Product Name	CD69 Rabbit mAb
Remarks/Alias	AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26; CD69
Species	Human
GeneID (Human)	969
GeneID	969
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 62-199 of human CD69 (NP_001772.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q07108
Protein Weight	23kDa
Shipping	Ice bag