

Anti-Cd74 Monoclonal Antibody For Wb, Ihc-P, Elisa - P04233

Item Number: CM0008983



Introduction

Rabbit monoclonal antibody against human CD74 (invariant chain) validated for WB, IHC-P, and ELISA. Targets a key MHC class II chaperone and MIF receptor, supporting immune research and antiviral studies.

[Learn More](#)

Parameter	Value
Product Name	CD74 Rabbit mAb
Remarks/Alias	II; p33; CLIP; DHLAG; HLADG; Ia-GAMMA; CD74
Species	Human
Gene ID (Human)	972
Gene ID	972
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD74 (P04233).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P04233
Protein Weight	34kDa
Shipping	Ice bag