

Anti-Cd79A Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P11912

Item Number: CM0006520



Introduction

Rabbit monoclonal antibody targeting CD79a (Ig-alpha), validated for WB, IHC, IF, and ELISA in human, mouse, and rat. Essential for B-cell antigen receptor signaling and B-cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	CD79a Rabbit PolymAb®
Remarks/Alias	IGA, MB1, MB-1, IGA1alpha
Species	Human
Gene ID (Human)	973
Gene ID	973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-226 of human CD79a. (NP_001774.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11912
Protein Weight	25kDa
Shipping	Ice bag