

Anti-Cd79B Monoclonal Antibody For Wb, Ihc-P, Elisa - P40259

Item Number: CM0002803



Introduction

Anti-CD79B rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Detects human and mouse CD79B (Ig-beta, B29), a B-cell receptor signaling protein essential for immunity. Suitable for B-cell studies.

[Learn More](#)

Parameter	Value
Product Name	CD79B Rabbit mAb
Remarks/Alias	B29; IGB; AGM6; Igbeta; CD79B
Species	Human
Gene ID (Human)	974
Gene ID	974
Immunogen	Synthetic peptide. A synthetic peptide corresponding to a sequence within amino acids 130-229 of human CD79B (NP_000617.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40259
Protein Weight	26kDa
Shipping	Ice bag