


Anti-Cd86 (B7-2) Polyclonal Antibody For Wb, If/Icc, Elisa - P42081

Item Number: CM0013827



Introduction

Rabbit polyclonal antibody against human CD86 (B7-2), validated for western blotting, immunofluorescence, and ELISA. Suitable for T-cell activation and immunology research applications.

[Learn More](#)

Parameter	Specification
Product Name	CD86 Rabbit pAb
Remarks/Alias	B70; B7-2; B7.2; LAB72; CD28LG2; CD86
Species	Human
GeneID (human)	942
GeneID	942
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-133 of human CD86 (NP_787058.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P42081
Protein Weight	38kDa
Shipping	Ice bag