


Anti-Cdc27 Monoclonal Antibody For Wb, Elisa - P30260

Item Number: CM0012764



Introduction

CDC27 rabbit monoclonal antibody validated for WB and ELISA, targeting human and mouse CDC27 (APC3). Key tool for cell cycle and APC/C ubiquitin ligase research, with a predicted weight of ~92kDa.

[Learn More](#)

Parameter	Value
Product Name	CDC27 Rabbit mAb
Remarks/Alias	APC3; HNUC; NUC2; H-NUC; ANAPC3; CDC27Hs; D0S1430E; D17S978E; CDC27
Species	Human
Gene ID (Human)	996
Gene ID	996
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 702-824 of human CDC27 (P30260).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30260
Protein Weight	92kDa
Shipping	Ice bag