

Anti-Cdc42Ep3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uki2

Item Number: CM0029467

Introduction

Rabbit polyclonal antibody against human CDC42EP3 (BORG2/CEP3), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Useful for studying actin cytoskeleton organization and cell morphology.

[Learn More](#)

Parameter	Value
Product Name	CDC42EP3 Rabbit pAb
Remarks/Alias	UB1; CEP3; BORG2; CDC42EP3
Species	Human
GeneID (Human)	10602
GeneID	10602
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-254 of human CDC42EP3 (NP_006440.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKI2
Protein Weight	28kDa
Shipping	Ice bag